

Assay Performance Characteristics
ulti med Products
www.ultimed.de

April 2004

Product Name: *DrugControl*[®] Barbiturates (BAR) 300 ng / mL

Product ID: Card 008A092

Dipstick 008T088

Cutoff Levels:

Assay	Calibrator	Cut-off (ng/mL)
Barbiturates (BAR)	Secobarbital	300

Specificity (sensitivity to chemically-related compounds):

	LEVEL (ng/mL)	% REACTIVITY
Barbiturates (BAR)		
Secobarbital	300	100%
Amobarbital	300	100%
Alphenol	150	200%
Aprobarbital	200	150%
Butabarbital	375	80%
Secobarbital	300	100%
Butalbital	2,500	12%
Butethal	100	300%
Cyclopentobarbital	600	50%
Pentobarbital	300	100%
Phenobarbital	100	300%