

Assay Performance Characteristics
ulti med Products
www.ultimed.de

April 2004

Product Name: *DrugControl*[®] Cocaine (COC) 150 ng / mL

Product ID: Card 008A116

Dipstick n / a

Cutoff Levels:

Assay	Calibrator	Cut-off (ng/mL)
Cocaine (COC)	Benzoyllecgonine	150

Specificity (sensitivity to chemically-related compounds):

	LEVEL (ng/mL)	% REACTIVITY
COCAINE (COC)		
Benzoyllecgonine	150	100%
Cocaine HCl	400	38%
Cocaethylene	6,250	2%
Ecgonine HCl	12,500	1%
Ecgonine methylester	50,000	<1%