

Assay Performance Characteristics  
ulti med Products  
www.ultimed.de

April 2004

Product Name: *DrugControl*<sup>®</sup> Metamphetamine (MET) 1000 ng / mL

Product ID: Card 008Z112

Dipstick 008T085

**Cutoff Levels:**

Assay	Calibrator	Cut-off (ng/mL)
Methamphetamine (MET)	D- Methamphetamine	1000

**Specificity (sensitivity to chemically-related compounds):**

	LEVEL (ng/mL)	% REACTIVITY
<b>METHAMPHETAMINE (MET)</b>		
$\rho$ -Hydroxymethamphetamine	30,000	3%
D-Methamphetamine	1,000	100%
L-Methamphetamine	8,000	12.5%
( $\pm$ )-3,4-Methylenedioxyamphetamine	2,000	50%
Mephentermine	50,000	2%