

Assay Performance Characteristics
ulti med Products
www.ultimed.de

April 2004

Product Name: *DrugControl*[®] Methadone (MTD) 300 ng / mL

Product ID: Card 008A090

Dipstick 008S100

Cutoff Levels:

Assay	Calibrator	Cut-off (ng/mL)
Methadone (MET)	Methadone	300

Specificity (sensitivity to chemically-related compounds):

	LEVEL (ng/mL)	% REACTIVITY
METHDONE (MTD)		
Methadone	300	100%
Doxylamine	50,000	<1%