

Assay Performance Characteristics
ulti med Products
www.ultimed.de

April 2004

Product Name: *DrugControl*[®] Phencyclidine (PCP) 25 ng / mL

Product ID: Card 008A093

Dipstick 008T087

Cutoff Levels:

Assay	Calibrator	Cut-off (ng/mL)
Phencyclidine (PCP)	Phencyclidine	25

Specificity (sensitivity to chemically-related compounds):

	LEVEL (ng/mL)	% REACTIVITY
Phencyclidine (PCP)		
4-Hydroxyphencyclidine	12,500	<1%
Phencyclidine	25	100%