

Assay Performance Characteristics  
ulti med Products  
www.ultimed.de

April 2004

Product Name: *DrugControl*<sup>®</sup> Marijuana (THC) 50 ng / mL

Product ID: Card 008Z115

Dipstick 008T083

**Cutoff Levels:**

Assay	Calibrator	Cut-off (ng/mL)
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9 COOH	50

**Specificity (sensitivity to chemically-related compounds):**

	LEVEL (ng/mL)	% REACTIVITY
<b>Marijuana (THC)</b>		
Cannabinol	20,000	<1%
11-nor- $\Delta^8$ -THC-9 COOH	30	167%
11-nor- $\Delta^9$ -THC-9 COOH	50	100%
$\Delta^8$ -THC	15,000	<1%
$\Delta^9$ -THC	15,000	<1%