

SalivaScreen®

EU-Patent • Not for sale in the US

Oral-Fluid-DOA Oral Fluid Drugtest

NEW FEATURE!

This saliva indicator** turns red, if the collection sponge is sufficiently saturated with saliva, to perform a test.


No more spoiled tests due to insufficient sample material!

**patent pending

Possible Drugparameters and Cut-Offs

AMP 50 ng/ml	EDDP 20 ng/ml	OXY 40 ng/ml
BUP 5 ng/ml	MET 50 ng/ml	PCP 10 ng/ml
BZD 10 ng/ml	MTD 30 ng/ml	THC 50* ng/ml
COC 20 ng/ml	OPI 40 ng/ml	THC 12 ng/ml

* THC Parent - Δ -9-tetrahydrocannabinol
Also available with Alcohol Teststrip at: ALC 0.02 %



①




②



③



④




⑤



GECKO Pharma proudly presents the new **SalivaScreen®** drugtesting device „**Oral-Fluid-DOA**“. This new saliva drugtest device comes with the following new features.



fig. A/B



- Saliva indicator
- Lock feature - the swab locks into the device for 100% tamper proof, airtight seal
- Extraction ports for lab confirmation (see figure A/B)
- Market leading drug cutoffs - No buffer required
- Detect up to 10 drugs + Alcohol
- One drug per strip - superior accuracy and sensitivity
- Crystal-clear container with magnification of result windows
- Advance sponge material
 - no additives
 - best recovery of THC
 - optimized for faster running tests

GECKO
PHARMA

GECKO Pharma GmbH • Reeshoop 1 • 22926 Ahrensburg • GERMANY

www.geckopharma.de • info@geckopharma.de