



Drug Test Cassette
008S901 AMP-OPI-THC/MTD-MET-COC

* protection by utility patents

Advantages

- reliable detection in saliva
- up to 6 drugs in just one stage
- saliva samples almost impossible to manipulate
- quick and easy
- immediate positive results on THC consumption
- black housing facilitates* interpretation through greater contrast to test field
- more distinct lines on negative results
- results will be stored on internal memory* (only with SalivaScreen Handheld Reader Touch)

Our service for you

Long time stability

In general, all ulti med products are stable for a minimum of 18 months after production.

Tests can be stored at room temperature, up to + 30°C.

High reliability

The ulti med rapid tests are extreme reliable and match in up to 99% (depending on test) with results of the reference method. They are also highly specific.

Competent advice

Our team will be at your disposal for any kind of question concerning the handling of our products.

Special questions concerning interaction and cross-reactivity with prescribed medicine cannot always be answered on the phone despite of our many years of experience. In this case please send your questions by fax or e-mail to our specialists. Usually you will receive a reply within 2 business days.

Quality

ulti med does not only fulfill the compulsory directions of their Quality Assurance System according to ISO EN 13485:2003, but also meets requirements which exceed these strict standards in particular with regard to extensive testing and complaint management.

Further information

You will find all information concerning our products on our website www.ultimed.org - of course, upon your request we will also send you any information by e-mail, fax or common mail.

Requests for product information please contact: info@ultimed.org

Research and Development / Production and Sales

Phone: +49 (0)4102-800 90

Fax: +49 (0)4102-500 82

e-mail: info@ultimed.org

Adress: ulti med Products (Deutschland) GmbH

Reeshoop 1

22926 Ahrensburg, Germany

P.O. Box 1814, 22908 Ahrensburg

Quality assurance for your safety

ulti med Products (Deutschland) GmbH complies with all requirements for production and sales of IVD (in vitro diagnostics) and is certified according to

EN ISO 13485:2003.

Tax ID-Nr. DE 154 23 11 80

Register of Companies / Amtsgericht Ahrensburg: HRB 3562

Notified body TÜV SÜD Munich Germany



SalivaScreen® Handheld Reader Touch
*Digital saliva drug test analyzer,
 with integrated Printer and
 3.5" high-resolution touch-screen display*



Saliva Collector

Specimen Collection

The oral fluid specimen should be collected using the collection device provided in the kit.

Follow the detailed directions for use.

No other collection devices should be used with this assay.

Oral fluid collected at any time of the day may be used.

The **SalivaScreen**[®]

Handheld Reader Touch

is a compact device that permits the digital interpretation and data management of saliva drug tests.

Therefore, it is the ideal device for roadside testing, emergency unit use, for emergency and corporate physicians, detention centres, border police etc.

- easy to use
- secure
- reliable
- complete with operating system
- optimized menu navigation with touch screen
- integrated Printer

Benefits:

- 3.5" touch-screen display
- Optimized touch panel navigation
- Detection of up to 6 drugs in human saliva in one test
- Direct and fast evaluation
- Integrated backup function. Evaluation results will optional be stored in device and testcassette as well.
- GPS position tracking system
- Fingerprint scanner
- Barcode reader
- 2 USB connections
- GPRS internal (optional)
- SD card reader (max. 32 GB)
- Card reader for magnetic strips and micro chips (i.e. Credit Cards)
- 4 x Sim Card Slots
- RFID – Radio Frequency Identification Device / tag STD



Prod. No.	Product
008M900	SalivaScreen Handheld Reader Touch®, incl. Alu-Transportation-Case and Leather-Bag <i>technically approved & certified for: EN 61326 and EN 61010</i>
008S901	SalivaScreen Drug Test Cassette (with built in memorychip)
008F120	Thermo-Print Paper roll

SalivaScreen[®]

- Drug detection test for saliva samples.


Detection of up to 6 drugs in one test / with internal memory to store the results of the evaluation

The following drugs can be detected:

- | | |
|--------------------------|-----------|
| ■ Amphetamines/MDA | 50 ng/ml |
| ■ Cocaine (Parent) | 20 ng/ml |
| ■ Methamphetamines/MDMA | 50 ng/ml |
| ■ Methadone | 30 ng/ml |
| ■ Marijuana (THC Parent) | 100 ng/ml |
| ■ Opiates | 40 ng/ml |

Area of application

The SalivaScreen Cassette is intended for use by professional institutions, drug clinics, police, border police, company employee screening programs and other point of contact test sites.



It is a useful tool for all operators as saliva samples provide certain advantages over other liquids like urine and blood. Drugs can be detected in saliva samples in a very reliable manner.

Saliva samples are easily accessible and can hardly be manipulated. Drugs can be detected immediately after consumption even before they appear in a metabolized state in urine. Basically, the SalivaScreen Cassette has been designed with lower cut-offs as drug tests designed for the examination of urine samples. Saliva samples show a far better correlation with blood samples than urine samples.