


ulti med Streptococcus A Test Cassette - Comparison of Flocked Swabs to Polyester Standard Swabs

The new swab in the ulti med Streptococcus A Card offers several advantages:

- Perpendicular Nylon fibers act like a soft brush and allow improved collection of cell samples.
- Capillary action between the Nylon fiber strands facilitates uptake of liquid samples.
- Sample stays close to the surface allowing easier elution.

Thus, the flocked swab increases the sensitivity of the Streptococcus A antigen test kit.

Detection of Streptococcus A Antigen*



* The comparison of the swabs was performed with the Streptococcus A Test Cassette and different dilutions of Streptococcus A antigen